

Prozomix Limited, Station Court, Haltwhistle, Northumberland, NE49 9HN, UK. Tel: +44 (0) 1434 400455 Fax: +44 (0) 1434 322822 info@prozomix.com technical@prozomix.com

[full product listing at www.prozomix.com/products/listing]

add this product to cart

view all branching enzymes

Recombinant Enzyme Product Specification Sheet

Cat. No.: PRO-E0406

www.prozomix.com

LOT: 2009-0406

Activity: 1,4-α-Glucan branching enzyme

Synonyms: Branching enzyme, amylo- $(1,4\rightarrow1,6)$ -transglycosylase; Q-enzyme; α -

glucan-branching glycosyltransferase; amylose isomerase;

enzymatic branching factor; branching glycosyltransferase; enzyme Q; glucosan transglycosylase; glycogen branching enzyme; plant branching enzyme; α -1,4-glucan: α -1,4-glucan-6-glycosyltransferase; starch branching enzyme; 1,4- α -D-glucan:1,4- α -D-glucan 6- α -D-(1,4- α -D-glucan)

α-D-glucano)-transferase

Nomenclature: CAZy [GH13 subf9, glycoside hydrolase family 13 subfamily 9,

member of clan GH-H]

Source organism: Escherichia coli str. K-12 substr. W3110

Enzyme Commission No.: 2.4.1.18

Activity: 45.58 U/ml

Specific activity: (37°C; pH 7.5; 3.3 mg/mL starch)

Purity: > 95 % as judged by SDS-PAGE

Form and storage: Supplied in 3.2 M ammonium sulphate, store at 4°C (shipped at room

temperature)

pH optimum: -

Temperature optimum: -

[Protein]: 2.952 mg/mL

Sequence length: 728 amino acids (view sequence)

Accession No.: BAE77860

Molecular weight: 88157.0 Da (theoretical)

~ 90000 Da (observed by SDS-PAGE)

(observed by mass spectrometry)

Biological function: Transfers a segment of a $(1->4)-\alpha$ -D-glucan chain to a primary

hydroxy group in a similar glucan chain

Potential application(s): Carbohydrate research, fundamental research



Prozomix Limited, Station Court, Haltwhistle, Northumberland, NE49 9HN, UK. Tel: +44 (0) 1434 400455 Fax: +44 (0) 1434 322822 info@prozomix.com technical@prozomix.com

www.prozomix.com

[full product listing at www.prozomix.com/products/listing]

Comments: Converts amylose into amylopectin. The description (official name)

requires a qualification depending on the product, glycogen or amylopectin, e.g. glycogen branching enzyme, amylopectin

branching enzyme. The latter has frequently been termed Q-enzyme

Usage: Agitate bottle sufficiently to fully homogenise enzyme precipitate

before use. Dilute in 50 mM sodium phosphate buffer, pH 7.5,

containing 2 mg/mL BSA. Do not dilute in water

Assay: One unit is defined as the amount of enzyme required to cause a fall

of 1.0 absorbance units, where the reaction mixture comprises 3.33 mg/mL starch (Sigma; S-9765) in 41.7 mM sodium phosphate buffer, pH 7.5, and where 0.24 mL of the reaction mixture is withdrawn at each time point, and mixed with 1.0 mL of deionised water and 0.2

mL of iodine reagent immediately prior to reading at 660 nm

Primary sequence:

MSDRIDRDVINALIAGHFADPFSVLGMHKTTAGLEVRALLPDATDVWVIEPKTGRKLAKLECLDSRGFFSGVIPR RKNFFRYQLAVVWHGQQNLIDDPYRFGPLIQEMDAWLLSEGTHLRPYETLGAHADTMDGVTGTRFSVWAPNARRV SVVGQFNYWDGRRHPMRLRKESGIWELFIPGAHNGQLYKYEMIDANGNLRLKSDPYAFEAQMRPETASLICGLPE KVVQTEERKKANQFDAPISIYEVHLGSWRRHTDNNFWLSYRELADQLVPYAKWMGFTHLELLPINEHPFDGSWGY QPTGLYAPTRRFGTRDDFRYFIDAAHAAGLNVILDWVPGHFPTDDFALAEFDGTNLYEHSDPREGYHQDWNTLIY NYGRREVSNFLVGNALYWIERFGIDALRVDAVASMIYRDYSRKEGEWIPNEFGGRENLEAIEFLRNTNRILGEQV SGAVTMAEESTDFPGVSRPQDMGGLGFWYKWNLGWMHDTLDYMKLDPVYRQYHHDKLTFGILYNYTENFVLPLSH DEVVHGKKSILDRMPGDAWQKFANLRAYYGWMWAFPGKKLLFMGNEFAQGREWNHDASLDWHLLEGGDNWHHGVQ RLVRDLNLTYRHHKAMHELDFDPYGFEWLVVDDKERSVLIFVRRDKEGNEIIVASNFTPVPRHDYRFGINQPGKW REILNTDSMHYHGSNAGNGGTVHSDEIASHGRQHSLSLTLPPLATIWLVREAE

Literature: 1. Hayashi et al. (2006) Mol. Syst. Biol. 2, 1-5